



# AFM Controlled Scanning Electrochemical Microscopy (AFM/SECM) With Seminal Nanonics Advances

Nanonics Imaging Ltd


[www.nanonics.co.il](http://www.nanonics.co.il)

# Straight Glass Pipette Probes with a Platinum Core Are Generally Accepted As The Best Probe for Scanning Electrochemical Microscopy (SECM)

***In spite of the Excellence of Glass Probes for SECM Several Issues Have Retarded Advances in SECM.***

***Principally these issues are:***

- The general straight geometry of the glass probes lack atomic force feedback capability
- Thus, the measured electrochemical current has to be used as a feedback mechanism
- As a result the measured current is compromised in its functional relevance due to alterations in current due to feedback
- Also the electrochemical current is not a sensitive feedback as is the case for atomic force microscopy
- Therefore, probe diameters in the 1-3 micron range are standard
- **Thus the excellence of glass probes is today seriously compromised in terms of SECM and the resolutions it can achieve**




# Nanonics' Has Transformed SECM By Resolving The Main Issue Impeding SECM Progress While Maintaining the Proven Excellence of Glass SECM Probes


## ***New Horizons in SECM With Nanonics.***

- Nanonics has 25 years of expertise in supplying SPM systems with exclusive glass ultrasensitive SPM probe technology
- Using this expertise Nanonics has transformed the excellence of glass SECM straight probes with Nanonics Systems into ultrasensitive AFM with SECM


Free optical axis from  
above and below




Single Probe  
Series



Hydra Multiprobe Series for  
Investigating Electron Transport  
with 2 SECM Probes





Such transparent integration with upright microscopes has been a hallmark of Nanonics systems for its establishment 25 years ago




Cantilevered Glass Based AFM Probes  
Exclusive for Nanonics' Systems

*What customers publish: " ...In our experimental configuration, the tip apex is not shadowed, as the laser beam is a focused beam with a converging angle of  $26.7^\circ$ , while the half angle of the tip is only about  $4^\circ$ . " [Sun and Shen, Apertureless near-field scanning Raman microscopy using reflection scattering geometry, Ultramicroscopy 94, 237 (2003)]*

**Compatibility With All  
Third Party AFM Probes**  
*But Such Probes In  
General Block The Optical  
Axis*



# And The Hydra MultiProbe Permits Probes To Be Brought Into Contact & Separated With Nanometric Precision *For Investigation of Nanotransport of Electron Transport, Conductivity, Optical & Thermal Phenomena*



**What customers publish:** " Collisions between both tips were effectively prevented by a proximity detection mechanism based on mechanical interaction between both tips. The minimum distance between both tips is reached when both tip coatings touch.." Klein et al, NanoLett. 14, 5010 (2014)




# Nanonics' Has Transformed SECM By Resolving The Main Issue Impeding SECM Progress While Maintaining the Proven Excellence of Glass SECM Probes

## ***New Horizons in SECM With Nanonics.***

- Nanonics has 25 years of expertise in supplying SPM systems with exclusive glass ultrasensitive SPM probe technology
- Using this expertise Nanonics has transformed the excellence of glass SECM straight probes with Nanonics Systems into ultrasensitive AFM with SECM
- Neither the systems nor the probes obstruct the optical axis from above and can be placed on any upright microscope

**Cantilevered Glass  
Based AFM Probes  
Exclusive for  
Nanonics' Systems**




Allowing The Probe To Be In Place With  
Optical Objectives with 100X 0.9NA  
1mm Working Distance  
and Water Immersion Objectives

# Nanonics' Has Transformed SECM By Resolving The Main Issue Impeding SECM Progress While Maintaining the Proven Excellence of Glass SECM Probes

## ***New Horizons in SECM With Nanonics:***

- Nanonics has 25 years of expertise in supplying SPM systems with exclusive glass ultrasensitive SPM probe technology
- Using this expertise Nanonics has transformed the excellence of glass SECM straight probes with Nanonics Systems into ultrasensitive AFM with SECM
- Neither the systems nor the probes obstruct the optical axis from above and can be placed on any upright microscope
- **The tuning fork AFM feedback methodology is used allowing for water immersion objectives and permits ultra high force sensitivity of the AFM with SECM down to pico Newtons**


AFM Sensing  
Cantilevered Glass  
SECM Probe with  
Tuning Fork






## **What customers publish:**

" measurements were carried out in non-contact mode at a constant separation distance of 4–8 nm from the surface using a Nanonics Multiview 2000 AFM equipped with a long AFM probe specially made for electrochemical measurements. The accurate separation distance is determined based on the error caused during tip positioning by a system generated feedback. The AFM probe tip with a frequency of 39.95 kHz and a Q factor (quality factor) of 1500 helps in providing better stability for performing scans in liquids." **Nanoscale, 2018, 10, 6962**

# 3D electrochemical deposition of Copper pillars with no jump to contact high rigidity quartz crystal tuning fork (QTF) nanopipettes



100:1  
Aspect  
Ratio



**Table 1.** The effects of the micropipette inner diameter on the printed Cu pillar diameter, deposition current, deposition current density, and average deposition rate.

Pipette size [μm]	Pillar diameter [μm]	Current, <i>I</i> [nA]	Current density [mA cm <sup>-2</sup> ]	Average deposition rate [nm s <sup>-1</sup> ]
1	1.5	5	283	90
1	1.5	6	340	105
0.4	0.55	1	421	160
0.4	0.61	1.2	411	130
0.1	0.23	0.3	722	300
0.1	0.27	0.4	699	230

## What customers publish:


“The development of a 3D electrochemical deposition system, which combines meniscus-confined electrodeposition (MCED) with atomic force microscope (AFM) closed-loop control and has a submicron resolution, is described. Thanks to the high rigidity of the hollow borosilicate glass (or quartz) tip and quartz crystal tuning fork (QTF), combined with the QTF’s high force sensitivity, the use of a solution-filled AFM tip in air is successful.” *Adv. Mater. Technol.* **2020**, 1900827



# Full Optical Integration Permits On-line Spectral Characterization Including Raman Chemical Identification With Electrochemical Analysis

## ***New Horizons in SECM With Nanonics:***


- Nanonics has 25 years of expertise in supplying SPM systems with exclusive glass ultrasensitive SPM probe technology
- Using this expertise Nanonics has transformed the excellence of glass SECM straight probes with Nanonics Systems into ultrasensitive AFM with SECM
- Neither the systems nor the probes obstruct the optical axis from above and can be placed on any upright microscope
- The tuning fork AFM feedback methodology is used allowing for water immersion objectives and permits ultra high force sensitivity of the AFM with SECM down to pico Newtons
- **This allows Nanonics to transparently unify SECM and Raman for electrochemical correlation with chemical analysis**



Raman  
MicroSpectroscopy

SECM  
Atomic Force Microscopy

Laser Excitation and  
Raman Scattering  
with SECM







AFM Raman

**The  
Raman**

**The  
Raman**


**The  
AFM**





As with all exclusive Nanonics probes the geometry of the probes and the instruments permit such unique applications in liquid like TERS and SECM with Raman

Laser Excitation and  
Raman Scattering  
with SECM



The above data was obtained with the probe inserted into a water immersion objective which is critical since without such an ability the Raman signal is low to non-existent

## *Come up to a new world of stand alone SECM by adding AFM*

- Get all standard modes of SECM such as SVET (Scanning vibrating electrode technique) and LEIS (localized electrochemical impedance spectroscopy) etc
- *But* with on-line AFM one has the ability to get not only micrometric but also nanometric imaging since AFM provides an independent feedback mode conventional SECM does not provide. This is essential to separate topography from electrochemical current which one needs to know independently.
- Thus, ultrasmall electrodes down to 100nm can be used which were previously impossible to employ
- Also allow important upgrades such as on-line Raman for chemical characterization along with ionic currents unavailable with previous SECM approaches
- Today SECM systems without AFM are missing the full picture without an independent feedback separating topography from electrochemical current and without accurate control of distances from the surface.